

EPO ELISA

FDA Cleared

NEW 0.6 mIU/mL LOW END SENSITIVITY

ERYTHROPOIETIN

Quantitative Measurement of Erythropoietin (EPO)

EPO: For the Diagnosis of Anemias and Erythrocytosis

Unique

Diagnostic adjunct in determining cause of anemia or erythrocytosis and Monitoring of the drug Rx Epogen.

Advanced

Direct label 2-site ELISA for reliability and convenience, capture antibody with distinctly different epitope is conjugated with biotin
Enhanced 0.6 mIU/mL low-end sensitivity

Accurate

Excellent patient sample correlation

Simple

Single step incubation and one washing step (2-site sandwich assay)

Convenient

Same-day Assay - 2 hour incubation for the sandwich assay+1/2 hour reaction color development



BIOMERICA

INTENDED USE

The Biomerica EPO Immunoassay is a two-site ELISA [Enzyme-Linked ImmunoSorbent Assay] for the quantitative measurement of the biologically active 165 amino acid chain of Erythropoietin (EPO) in human serum. The simple 2 ½ hour ELISA test is intended for in vitro diagnostic use, as an aid in the diagnosis of anemias and polycythemias. With the advent of the administration of recombinant erythropoietin as a biologic therapy to increase red blood cell mass, an erythropoietin assay may be used also to aid in the prediction and monitoring of response to recombinant erythropoietin treatment in persons with anemias.

BACKGROUND

Human EPO is a polypeptide consisting of 165 amino acids and is a heavily glycosylated protein with a molecular weight of about 30,000 - 34,000 Daltons. Serum EPO levels are dependent on the rate of production and the rate of clearance of the protein. Ninety percent of EPO is produced in the peritubular cells of the adult kidney in response to a decrease in tissue oxygenation^{3,4}.

Quantitation of serum erythropoietin concentration serves as a diagnostic adjunct in determining the cause of anemia or erythrocytosis. Aplastic anemia, hemolytic anemia and anemia due to iron deficiency all result in serum EPO elevation. Whereas, EPO levels in patients with secondary anemia due to renal failure and other disorders such as acquired immune deficiency syndrome (AIDS) are generally inappropriately low for the degree of anemia. The illegal use of EPO and EPO derivatives by athletes has also spurred demand for EPO testing.

Polycythemia rubra vera, or primary erythrocytosis (an increase of red blood cell mass) results from unstimulated over production of erythrocytes. Hence, the increase in the hemoglobin causes decreased production of EPO, which results in subnormal levels of serum EPO⁹.

PERFORMANCE

	<u>Assay Time</u>	<u>Accuracy</u>	<u>Sensitivity</u>
EPO:	2 1/2 hrs.	Biomerica ELISA = 0.94 ELISA Kit - 0.41 mIU/mL r=0.989 N=85	0.6 mIU/mL

ORDERING

Catalog No.	Description
7025	EPO ELISA kit - Quantitative (96 tests)
7025B	EPO Bulk - Unlabeled (minimum quantities required)

CE and EN ISO 13485:2003 Compliant, Multi-language inserts available

Bibliography

1. Sawyer, S.T., Krantz, S.B., Sawada, K. Receptors for Erythropoietin in Mouse and Human Erythroid Cells and Placenta. *Blood* 1989; 74: 103-109.
2. Imai, N., Kawamura, A., Higuchi, M., et al. Physicochemical and Biological Comparison of Recombinant Human Erythropoietin with Human Urinary Erythropoietin. *J Biochem* 1990; 107: 352-359.
3. Jacobson, L.O., Goldwasser, E., Fried, W., Pizak, L.F. The Role of the Kidney in Erythropoiesis. *Nature* 1957; 179: 633-634.
4. Koury, S.T., Bondurant, M.C., Koury, M.J. Localization of Erythropoietin Synthesizing Cells in Murine Kidney by in-situ Hybridization. *Blood* 1988; 71: 524-527.
5. Goldberg, M.A., Dunning, S.P., Bunn, H.F. Regulation of the Erythropoietin Gene: Evidence that the Oxygen Sensor is a Heme Protein. *Science* 1988; 242: 1412-1415.
6. Erslev, A.J., Caro, J., Birgegard, G., Silver, R., Miller, O. The Biogenesis of Erythropoietin. *Experimental Hematology* 1980; Suppl 8: 1-13.
7. Spivak, J.L. The Mechanism of Action of Erythropoietin. *Int J Cell Cloning* 1986; 4: 139-166.
8. Erslev, A.J. Erythropoietin. *New Eng J Med* 1991; 324:1339-1344.
9. Garcia, J.F., Ebbbe, S.N., Hollander, L., Cutting, H.O., Miller, M.E., Cronkites, E.P. Radioimmunoassay of Erythropoietin: Circulating Levels in Normal and Polycythemic Human Beings. *J Lab Clin Med* 1982; 99: 624-635.
10. Wild, D., editor. *The Immunoassay Handbook*, Stockton Press, 1994, p. 428.
11. Wide L., Bengtsson C., Birgegard G. Circadian Rhythm of Erythropoietin in Human Serum. *Br J Haematol* 1989; 72: 85-90.
12. Cahan C., Decker M.J., Arnold J.L., Washington L.H., Veldhuis J.D., Goldwasser E., Strohl K.P. Diurnal Variations in Serum Erythropoietin Levels in Healthy Subjects and Sleep Apnea Patients. *J Appl Physiol* 1992; 72: 2112-7.
13. Goldwasser E. and Sherwood J.B. Annotation, Radioimmunoassay of Erythropoietin. *Br J Haematol* 1981; 48: 359-63.
14. Kricka L.J. Human Anti-Animal Antibody Interferences in Immunological Assays. *Clin Chem* 1999; 45: 942-956.



BIOMERICA

17571 Von Karman Avenue • Irvine, CA 92614 USA • Tel (949) 645-2111
Fax (949) 553-1231 • email: bmra@biomerica.com • www.biomerica.com