

One Step Barbiturates Test Device (Urine) Package Insert

(Catalog Number: 1184-C)



A rapid, one step test for the qualitative detection of Barbiturates in human urine.

For in vitro diagnostic use only.

INTENDED USE

The One Step Barbiturates Test Device (Urine) is a lateral flow chromatographic immunoassay for the detection of Barbiturates (Secobarbital) in urine at a cut-off concentration of 300 ng/mL.

This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

SUMMARY

Barbiturates are central nervous system depressants. They are used therapeutically as sedatives, hypnotics, and anticonvulsants. Barbiturates are almost always taken orally as capsules or tablets. The effects resemble those of intoxication with alcohol. Chronic use of Barbiturates leads to tolerance and physical dependence. Short acting Barbiturates taken at 400 mg/day for 2-3 months produces a clinically significant degree of physical dependence. Withdrawal symptoms experienced during periods of drug abstinence can be severe enough to cause death. Only a small amount (less than 5%) of most Barbiturates are excreted unaltered in the urine. The detection period for the Barbiturates in the urine is 4-7 days.

The One Step Barbiturates Test Device (Urine) is a rapid urine-screening test that can be performed without the use of an instrument. The test utilizes a monoclonal antibody to selectively detect elevated levels of Barbiturates in urine. The One Step Barbiturates Test Device (Urine) yields a positive result when the Barbiturates in urine exceeds cut-off level.

PRINCIPLE

The One Step Barbiturates Test Device (Urine) is an immunoassay based on the principle of competitive binding. Drugs that may be present in the urine specimen compete against the drug conjugate for binding sites on the antibody.

During testing, a urine specimen migrates upward by capillary action. Barbiturates, if present in the urine specimen below the detectable cut-off level, will not saturate the binding sites of the antibody in the test. The antibody coated particles will then be captured by immobilized Barbiturates-protein conjugate and a

visible colored line will show up in the test line region. The colored line will not form in the test line region if the Barbiturates level exceeds the detectable cut-off level because it will saturate all the binding sites of anti-Barbiturates antibody.

A drug-positive urine specimen will not generate a colored line in the test line region because of drug competition, while a drug-negative urine specimen will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

REAGENTS

The test device contains anti-Barbiturates antibody coupled particles and Barbiturates-protein conjugate. A goat antibody is employed in the control line system.

PRECAUTIONS

- *For in vitro diagnostic use only.* Do not use after the expiration date.
- The test device should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The test device should be discarded according to federal, state and local regulations.

STORAGE AND STABILITY

The kit can be stored at room temperature or refrigerated (2-30°C). The test device is stable through the expiration date printed on the sealed pouch. The test device must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

SPECIMEN COLLECTION AND PREPARATION

Urine Assay

The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible particles should be centrifuged, filtered, or allowed to settle to obtain a clear specimen for testing.

Specimen Storage

Urine specimens may be stored at 2-8°C for up to 48 hours prior to testing. For long-term storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.

MATERIALS

Materials Provided

- Test devices
- Disposable droppers
- Package insert


Materials Required But Not Provided

- Specimen collection container
- Timer

DIRECTIONS FOR USE

Allow the test device, urine specimen and/or controls to reach room temperature (15-30°C) prior to testing.

1. Bring the pouch to room temperature before opening it. Remove the test device from the sealed pouch and use it as soon as possible.
2. Place the test device on a clean and level surface. Hold the dropper vertically and **transfer 3 full drops of urine** (approx. 100µl) to the specimen well (S) of the test device, and then start the timer. Avoid trapping air bubbles in the specimen well (S). See the illustration below.
3. Wait for the red line(s) to appear. The result should be **read at 5 minutes.** Do not interpret the result after 10 minutes.



INTERPRETATION OF RESULTS

(Please refer to illustration above)

NEGATIVE:* **Two lines appear.** One red line should be in the control region (C), and another apparent red or pink line should be in the test region (T). This negative result indicates that the Barbiturates concentration is the detectable cut-off level.

* **NOTE:** The shade of red in the test region (T) may vary, but it should be considered negative whenever there is even a faint pink line.

POSITIVE: **One red line appears in the control region (C).** No line appears in the test region (T). This positive result indicates that the Barbiturates concentration exceeds the detectable cut-off level.

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test device. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

QUALITY CONTROL

A procedural control is included in the test. A red line appearing in the control region (C) is considered an internal positive procedural control. It confirms sufficient specimen

volume, adequate membrane wicking and correct procedural technique.

Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as good laboratory testing practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

1. The One Step Barbiturates Test Device (Urine) provides only a preliminary analytical result. A secondary analytical method must be used to obtain a confirmed result. Gas chromatography and mass spectrometry (GC/MS) are the preferred confirmatory methods.^{1,2,3}
2. The One Step Barbiturates Test Device (Urine) is a qualitative screening assay and cannot determine either the drug concentration in the urine or the level of intoxication.
3. It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
4. Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
5. A positive result indicates presence of the drug or its metabolites but does not indicate level or intoxication, administration route or concentration in urine.
6. A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.

PERFORMANCE CHARACTERISTICS

Accuracy

The accuracy of the One Step Barbiturates Test Device (Urine) was evaluated by testing clinical urine specimens. These clinical urine specimens were previously collected from subjects presenting for Drug Screen Testing. Presumptive positive results were confirmed by GC/MS. The following results were tabulated:

Method	GC/MS		Total Results
	Positive	Negative	
Biomerica One Step Test Device	Positive	3	45
	Negative	219	222
Total Results		45	267
% Agreement with GC/MS Analysis		93%	98%

Specificity

The following table lists compounds that are positively detected in urine by the One Step Barbiturates Test Device (Urine) at 5 minutes.

Compound	Concentration (ng/mL)
Secobarbital	300
Amobarbital	300

Aprobarbital	200
Butalbital	2500
Butabarbital	75
Phenobarbital	100
Pentobarbital	300
Cyclopentobarbital	600
Alphenol	150
Butethal	100

Effect of Urinary Specific Gravity

Fifteen (15) urine specimens of normal, high, and low specific gravity ranges were spiked with 150 ng/mL and 450 ng/mL of Secobarbital respectively. The One Step Barbiturates Test Device (Urine) was tested in duplicate using the fifteen neat and spiked urine specimens. The results demonstrate that varying ranges of urinary specific gravity does not affect the test results.

Effect of the Urinary pH

The pH of an aliquoted negative urine pool was adjusted to a pH range of 5 to 9 in 1 pH unit increments and spiked with Secobarbital to 150 ng/mL and 450 ng/mL. The spiked, pH-adjusted urine was tested with the One Step Barbiturates Test Device (Urine) in duplicate and interpreted according to the package insert. The results demonstrate that varying ranges of pH does not interfere with the performance of the test.

Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in urine not associated with Barbiturates. The following compounds show no cross-reactivity when tested with the One Step Barbiturates Test Device (Urine) at a concentration of 100 ng/mL.

Non Cross-Reacting Compounds

4-Acetamidophenol	Fenoprofen	Papaverine
Acetophenetidin	Furosemide	Penicillin-G
N-Acetylprocainamide	Gentisic acid	Pentobarbital
Acetylsalicylic acid	Hemoglobin	Perphenazine
Aminopyrine	Hydralazine	Phenelzine
Amitypyline	Hydrochlorothiazide	Phenobarbital
Amobarbital	Hydrocodone	Phentermine
Amoxicillin	Hydrocortisone	L-Phenylephrine
Ampicillin	O-Hydroxyhippuric acid	β-Phenylethylamine
Ascorbic acid	p-Hydroxy-methamphetamine	Phenylpropanolamine
D,L-Amphetamine	3-Hydroxytyramine	Prednisolone
Apomorphine	Ibuprofen	Prednisone
Aspartame	Atropine	Procaine
Benzilic acid	Iproniazid	Promazine
Benzoic acid	(±) - Isoproterenol	Promethazine
Benzphetamine	Isoxsuprine	D,L-Propranolol
Bilirubin	Ketamine	D-Propoxyphene
Brompheniramine	Ketoprofen	D-Pseudoephedrine
Caffeine	Labetalol	Quinidine
Cannabinol	Levorphanol	Quinine
Chloralhydrate	Loperamide	Ranitidine
Chloramphenicol	Maprotiline	Salicylic acid
		Serotonin (5-

Chlordiazepoxide	Meperidine	Hydroxytyramine)
Chlorothiazide	Meprobamate	Sulfamethazine
(±) Chlorpheniramine	Methoxyphenamine	Sulindac
Chlorpromazine	(+) 3,4-Methylenedioxy-amphetamine	Temazepam
Chlorquine	(+) 3,4-Methylenedioxy-methamphetamine	Tetracycline
Cholesterol	Morphine-3-β-D glucuronide	Tetrahydrocortisone, 3 Acetate
Clomipramine	Morphine Sulfate	Tetrahydrocortisone 3 (β-D glucuronide)
Clonidine	Nalidixic acid	Tetrahydrozoline
Codeine	Naloxone	Thebaine
Cortisone	Naltrexone	Thiamine
(-) Cotinine	Naproxen	Thioridazine
Creatinine	Niacinamide	D, L-Tyrosine
Deoxycorticosterone	Nifedipine	Tolbutamide
Dextromethorphan	Norcodein	Triamterene
Diazepam	Norethindrone	Trifluoperazine
Diclofenac	D-Norpropoxyphene	Trimethoprim
Diflunisal	Doxylamine	Trimipramine
Digoxin	D,L-Octopamine	Tryptamine
Diphenhydramine	Oxalic acid	D, L-Tryptophan
Doxylamine	Oxalopam	Tyramine
Ecgonine methylester	Oxazepam	Uric acid
(-) Y Ephedrine	Oxalonic acid	Verapamil
Erythromycin	Oxycodone	Zomepirac
β-Estradiol	Oxymetazoline	
Estrone-3-sulfate		
Ethyl-p-aminobenzoate		

BIBLIOGRAPHY

1. Tietz NW. *Textbook of Clinical Chemistry*. W.B. Saunders Company, 1986; 1735
2. Baselt RC. *Disposition of Toxic Drugs and Chemicals in Man*. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488
3. Hawks RL, CN Chiang. *Urine Testing for Drugs of Abuse*. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986

	Storage Temperature		Authorized Representative
	Lot Code		Caution, See Instructions
	Expiration		For in vitro diagnostic use
	Manufacturer		Catalog No.



Biomerica, Inc., 1533 Monrovia Avenue
Newport Beach, CA 92663 USA



according to IVDD 98/79/ EC
MDSS
Burckhardtstrasse 1
30163 Hannover, Germany

