

One Step Methadone Test Device (Urine)

Package Insert

(Catalog Number: 1187-C)



A rapid, one step test for the qualitative detection of Methadone in human urine.

For *in vitro* diagnostic use only.

INTENDED USE

The One Step Methadone Test Device (Urine) is a lateral flow chromatographic immunoassay for the detection of Methadone in human urine.

This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

SUMMARY

Methadone is a narcotic pain reliever for medium to severe pain. It is also used in the treatment of Heroin (Opiate dependence: Vicodin, Percocet, Morphine, ect) addiction. Oral Methadone is very different than the IV Methadone. Oral Methadone is partially stored in the liver for later use. IV Methadone acts more like Heroin. In most states you must go to a pain clinic or a Methadone maintenance clinic to be prescribed Methadone. Methadone is a long acting pain reliever producing effects that last between twelve to forty-eight hours. Ideally, Methadone frees the client from the pressures of obtaining illegal Heroin, from the dangers of injection, and from the emotional roller coaster that most Opiates produce. Methadone if taken for long periods and at large doses can lead to a very long withdrawal period. As compared to other Opiates in which the withdrawal period is a week to ten days, heavy Methadone users can expect too not recover for up to 5 or 6 weeks.

The One Step Methadone Test Device (Urine) is a rapid urine-screening test that can be performed without the use of an instrument. The test utilizes a monoclonal antibody to selectively detect elevated levels of Methadone in urine. The One Step Methadone Test Device (Urine) yields a positive result when the Methadone in urine exceeds 300 ng/mL.

PRINCIPLE

The One Step Methadone Test Device (Urine) is an immunoassay based on the principle of competitive binding. Drugs that may be present in the urine specimen compete against the drug conjugate for binding sites on the antibody. During testing, a urine specimen migrates upward by capillary action. Methadone, if present in the urine specimen below 300 ng/mL, will not saturate the binding sites of antibody-coated particles in the test. The antibody coated particles will then be captured by immobilized Methadone-protein conjugate and a visible colored line will show up in the test line region. The

colored line will not form in the test line region if the Methadone level exceeds 300 ng/mL because it will saturate all the binding sites of anti-Methadone antibody. A drug-positive urine specimen will not generate a colored line in the test line region because of drug competition, while a drug-negative urine specimen will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

REAGENTS

The test device contains anti-Methadone antibody coupled particles and Methadone-protein conjugate. A goat antibody is employed in the control line system.

PRECAUTIONS

- For *in vitro* diagnostic use only. Do not use after the expiration date.
- The test device should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test device should be discarded according to federal, state and local regulations.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test device is stable through the expiration date printed on the sealed pouch. The test device must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

SPECIMEN COLLECTION AND PREPARATION

Urine Assay

The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitates should be centrifuged, filtered, or allowed settle to obtain a clear supernatant for testing.

Specimen Storage

Urine specimens may be stored at 2-8°C for up to 48 hours prior to assay. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.

PROCEDURE

Materials Provided

- Test devices
- Disposable specimen droppers
- Package insert

Materials Required But Not Provided

- Specimen collection container
- Timer


DIRECTIONS FOR USE

Allow test device, urine specimen, and/or controls to reach room temperature (15-30°C) prior to testing.

1. Bring the pouch to room temperature before opening it. Remove the test device from the sealed pouch and use it as soon as possible.
2. Place the test device on a clean and level surface. Hold the dropper vertically and **transfer 3 full drops of urine** (approx. 100µl) to the specimen well (S) of the test device,

and then start the timer. Avoid trapping air bubbles in the specimen well (S). See the illustration below.

3. Wait for the red line(s) to appear. The result should be **read at 5 minutes.** Do not interpret the result after 10 minutes.



INTERPRETATION OF RESULTS

(Please refer to the illustration above)

NEGATIVE:* Two lines appear. One red line should be in the control region (C), and another apparent red or pink line should be in the test region (T). This negative result indicates that the Methadone concentration is below the detectable cut-off level (300 ng/mL).

* **NOTE:** The shade of red in the test line region (T) will vary, but it should be considered negative whenever there is even a faint pink line.

POSITIVE: One red line appears in the control region (C). No line appears in the test region (T). This positive result indicates that the Methadone concentration exceeds the detectable cut-off level (300 ng/mL).

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test device. If the problem persists, discontinue using the lot immediately and contact your local distributor.

QUALITY CONTROL

A procedural control is included in the test. A red line appearing in the control region (C) is considered an internal positive procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as good laboratory testing practice to confirm the test procedure and to verify proper test performance.

LIMITATION

- The One Step Methadone Test Device (Urine) provides only a preliminary analytical result. A secondary analytical method must be used to obtain a confirmed result. Gas chromatography and mass spectrometry (GC/MS) are the preferred confirmatory methods.^{1,2}
- The One Step Methadone Test Device (Urine) is a qualitative screening assay and cannot determine either the drug concentration in the urine or the level of intoxication.
- It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
- Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
- A positive result indicates presence of the drug or its metabolites but does not indicate level or intoxication, administration route or concentration in urine.
- A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cutoff level of the test.

PERFORMANCE CHARACTERISTICS

Accuracy

The accuracy of One Step Methadone Test Device (Urine) was evaluated by testing clinical urine specimens. These clinical urine specimens were previously collected from subjects presenting for Drug Screen Testing. Presumptive positive results were confirmed by GC/MS. The following results were tabulated:

Method	Results	GC/MS		Total Results
		Positive	Negative	
Biomerica One Step Test Device	Positive	115	2	117
	Negative	2	181	183
	Total Results	117	183	300
% Agreement with GC/MS Analysis		98%	99%	98%

Analytical Sensitivity

A drug-free urine pool was spiked with (±) Methadone Hydrochloride at the following concentrations: 0 ng/mL, 150 ng/mL, 225 ng/mL, 300 ng/mL, 375 ng/mL and 450 ng/mL. The result demonstrates 100% accuracy at 50% above and 50% below the cut-off concentration. The data are summarized below:

(±) Methadone Hydrochloride Concentration (ng/mL)	Percent of Cut-off	n	Visual Result	
			Negative	Positive
0	0	30	30	0
150	-50%	30	30	0
225	-25%	30	26	4
300	Cut-off	30	18	12
375	+25%	30	5	25
450	+50%	30	0	30

Analytical Specificity

The following table lists compounds that are positively detected in urine by the One Step Methadone Test Device (Urine) at 5 minutes.

Compound

Concentration (ng/mL)

Methadone	300
Doxylamine	50,000

Effect of Urinary Specific Gravity

Fifteen (15) urine samples of normal, high, and low specific gravity ranges were spiked with 150 ng/ml and 450 ng/ml of Methadone respectively. The One Step Methadone Test Device (Urine) was tested in duplicate using the fifteen neat and spiked urine samples. The results demonstrate that varying ranges of urinary specific gravity does not affect the test results.

Effect of the Urinary pH

The pH of an aliquoted negative urine pool was adjusted to a pH range of 5 to 9 in 1 pH unit increments and spiked with Methadone to 150 ng/ml and 450 ng/ml. The spiked, pH-adjusted urine was tested with the One Step Methadone Test Device (Urine) in duplicate and interpreted according to the package insert. The results demonstrate that varying ranges of pH does not interfere with the performance of the test.

Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in urine not associated with Methadone. The following compounds show no cross-reactivity when tested with the One Step Methadone Test Device (Urine) at a concentration of 100 µg/mL.









Non Cross-Reacting Compounds


4-Acetamidophenol	Erythromycin	Papaverine
Acetophenetidin	β-Estradiol	Penicillin-G
N-Acetylprocainamide	Estrone-3-sulfate	Pentobarbital
Acetylsalicylic acid	Ethyl-p-aminobenzoate	Perphenazine
Aminopyrine	Flenfluramine	Phencyclidine
Amitypyline	Fenoprofen	Phenelzine
Amobarbital	Furosemide	Phenobarbital
Amoxicillin	Gentisic acid	Phentermine
Ampicillin	Hemoglobin	L-Phenylephrine
Ascorbic acid	Hydralazine	β-Phenylethylamine
D-Amphetamine	Hydrochlorothiazide	Phenylpropanolamine
D,L-Amphetamine	Hydrocodone	Prednisolone
L-Amphetamine	Hydrocortisone	Prednisone
Apomorphine	p-Hydroxyamphetamine	Procaine
Aspartame	O-Hydroxyhippuric acid	Promazine
Atropine	3-Hydroxytyramine	Promethazine
Benzilic acid	Ibuprofen	D,L-Propranolol
Benzoic acid	Imipramine	D-Propoxyphene
Benzoylcegonine	Iproniazid	D-Pseudoephedrine
Benzphetamine	(-) Isoproterenol	Quinacrine
Bilirubin	Isoxsuprine	Quinidine
Brompheniramine	Ketamine	Quinine
Caffeine	Ketoprofen	Ranitidine
Cannabidiol	Labetalol	Salicylic acid
Chloralhydrate	Levorphanol	Secobarbital
Chloramphenicol	Loperamide	Serotonin (5-Hydroxytyramine)
Chlordiazepoxide	Maprotiline	Sulfamethazine
Chlorothiazide	Meperidine	Sulindac
(±) Chlorpheniramine	Meprobamate	Temazepam
Chlorpromazine	Methoxyphenamine	Tetracycline
Chlorquine	(+) 3,4-Methylenedioxy-amphetamine	
Cholesterol		

Clomipramine	3,4-Methylenedioxyethyl-	Tetrahydrocortisone, 3
Clonidine	amphetamine	Acetate
Cocaethylene	Methylphenidate	Tetrahydrocortisone 3 (β-D glucuronide)
Cocaine hydrochloride		
Codeine	Morphine-3-β-D-glucuronide	Tetrahydrozoline
Cortisone		Thioridazine
(-) Cotinine	Nalidixic acid	D, L-Tyrosine
Creatinine	Naloxone	Tolbutamine
Deoxycorticosterone	Naltrexone	Trans-2-phenylcyclopropylamine
Dextromethorphan	Naproxen	Triamterene
Diazepam	Niacinamide	Trifluoperazine
Diclofenac	Nifedipine	Trimethoprim
Diffunilal	Norethindrone	Trimipramine
Digoxin	D-Norpropoxyphene	Triptamine
Diphenhydramine	Noscapine	D, L-Tryptophan
Doxylamine	D,L-Octopamine	Tyramine
Ecgonine hydrochloride	Oxalic acid	Uric acid
Ecgonine methylester (1R,2S)-(-)-Ephedrine	Oxazepam	Verapamil
L-Epinephrine	Oxolinic acid	Zomepirac
(-) Y Ephedrine	Oxycodone	
	Oxymetazoline	

BIBLIOGRAPHY

- Baselt RC. *Disposition of Toxic Drugs and Chemicals in Man*. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488
- Hawwks RL, CN Chiang. *Urine Testing for Drugs of Abuse*. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986

	Storage Temperature		Authorized Representative
	Lot Code		Caution, See Instructions
	Expiration		For in vitro diagnostic use
	Manufacturer		Catalog No.

 Biomerica, Inc., 1533 Monrovia Avenue
Newport Beach, CA 92663 USA



according to IVDD 98/79/ EC
MDSS
Burckhardtstrasse 1
30163 Hannover, Germany